





KERLIXTM AMDTM

ANTIMICROBIAL GAUZE

- **Low Cost** — Pennies per day
- **Convenient** — Works in moist or dry conditions
- **Simple** — Compliant with existing protocol



ZONE OF INHIBITION



- **KERLIX AMD Gauze vs. untreated gauze against Staphylococcus aureus.***
- **KERLIX AMD Gauze vs. untreated gauze against Escherichia coli.***
- **KERLIX AMD Gauze vs. untreated gauze against Candida albicans.***

* Dressing in contact with seeded agar plate modulated at 10^6 CFU/ml.

Kendall AMD Dressings




Kendall AMD Dressings can help control infection in two ways:

- **Primary Bacterial Barrier** —
AMD works as a contact kill and maintains log counts in the wound without negative impact to wound healing.
- **Secondary Layer** —
AMD works to limit cross-contamination by significantly reducing bacterial growth in the pooling layer.

AMD is effective against:

- Staphylococcus aureus
- Staphylococcus epidermidis
- Pseudomonas aeruginosa
- Escherichia coli
- Candida albicans
- Staphylococcus coagulase
- Proteus mirabilis
- Serratia marcescens
- Enterbacter cloacae
- Klebsiella pneumoniae
- Enterococcus faecalis

Graph: AMD Efficacy



*0/24/48 at 10⁶ inoculation levels for challenge organisms

Ordering Information

KERLIX AMD Super Sponges

Code	Description	Per Tray	Ship Case
Sterile 2's in Soft Pouch Packaging			
6662	Medium	40	480
Sterile 5's in Soft Pouch Packaging			
6665	Medium	50	600
Sterile 10's in Plastic Tray			
6660	Medium		480

KERLIX AMD Rolls

Code	Description	Ship Case
Sterile in Rigid Plastic Tray 6 ply, 1 roll/tray		
3331	4.5" x 4.1 yds	60 rolls
Sterile in Soft Pouch 6 ply, 1 roll/pkg		
3332	4.5" x 4.1 yds	100 rolls

™ Trademarks of Tyco Healthcare Group LP

tyco
Healthcare

Kendall

15 Hampshire Street
Mansfield, MA 02048
1-800-962-9888
508-261-8000

www.kendallhq.com

For International locations please see
www.tycohealthcare.com

H-5654 10M © Tyco 9/2004